



ZHIVAS LTD

INSTRUCTIONS FOR USE

ZHIVASEPT RAPID S

READY FOR USE ALCOHOL SOLUTION
FOR FAST DISINFECTION OF SMALL AND
HARDLY ACCESSIBLE SURFACES



Application

Zhivasept Rapid S is intended for rapid disinfection of alcohol-resistant surfaces and objects, including small and hard-to-reach surfaces, in health and medical institutions, in food production and trade enterprises, public facilities, in all industries and in households.

Properties

Zhivasept Rapid S is ready for use solution, aldehyde, phenol and quaternary ammonium compounds free. It contains a balanced synergic combination of alcohols, which guarantees fast and effective disinfection. Zhivasept Rapid S contains corrosion inhibitors, does not affect alcohol-resistant surfaces and is compatible with various materials.

The aerosol form of the preparation allows easy treatment of small and hard-to-reach surfaces.

Appearance Transparent liquid with pleasant aroma.
pH 5.5 – 7.0

Ingredients

Each 100 g ZHIVASEPT Rapid S contain

Active substances: 45 g Propan-1-ol (CAS 71-23-8)
 30 g Propan-2-ol (CAS 67-63-0)
 5 g Ethanol (CAS 64-17-5)

Auxiliary ingredients: Corrosion inhibitors, fragrance

Spectrum of activity

Bactericidal (incl. MRSA), mycobactericidal (incl. Tbc), virucidal (Adeno, Noro, Vaccinia – enveloped viruses) and fungicidal.

Microbiological activity	Test method	Exposition
Bactericidal (<i>P.aeruginosa</i> , <i>S.aureus</i> , <i>E.coli</i> , <i>E.hirae</i>)	EN 1276	1 min
Bactericidal (<i>P.aeruginosa</i> , <i>S.aureus</i> , <i>E.coli</i> , <i>E.hirae</i>)	EN 13697 / EN 13727	1 min
Bactericidal (MRSA)	EN 1040	
Fungicidal (<i>C.albicans</i> , <i>A.niger</i>)	EN 1650	1 min
Fungicidal (<i>C.albicans</i> , <i>A.niger</i>)	EN 13697 / EN13624	1 min
Mycobactericidal (<i>M.terrae</i> , <i>M.avium</i>)	EN 14348 / EN 14563	1 min
Virucidal (<i>Adeno</i> , <i>Noro</i> , <i>Vaccinia</i> - enveloped viruses)	EN 14476	1 min



ZHIVAS LTD

Method of application

Zhivasept Rapid S is applied undiluted. Disinfection is carried out by spraying the treated surfaces from a distance of about 30 cm, following the recommended exposition time shown below.

Disinfection object	Activity	Exposition
Various alcohol resistant surfaces in hospitals, healthcare institutions, public and commercial sector, food processing industry and households.	Bactericidal (<i>incl. MRSA</i>), mycobactericidal and fungicidal	1 minute
	Virucidal	1 minutes

Attention! Do not apply on or near appliances connected to the mains. Do not apply to surfaces near heat sources or open flames.

Treated surfaces and objects can be used immediately after the exposition time without subsequent rinsing with water.

Storage conditions

Keep out of the reach of children. Store in closed original packaging, in dry and well ventilated premises, away from sources of ignition, open flame heating devices, at $t < 30^{\circ} \text{C}$. Do not expose to direct sunlight.

Precautions

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If swallowed, do not induce vomiting, immediately seek medical advice and show this container or label.

Treatment of waste and packaging – according to local regulations

Do not release the solution in the environment – for disposal of waste, residuals and empty packaging, see the special instructions in the safety data sheet.

Packaging

Plastic bottles with sprayer pump of 750 ml and 1 L

Plastic canisters of 5 L.

Shelf life: 3 years

Manufacturer:

ZHIVAS Ltd, 1592 Sofia, Bulgaria
14, Assen Yordanov Blvd.
Tel/fax: +359 2 9817823, +359 2 9802507
zhivas@techno-link.com; www.zhivas.com