



INSTRUCTION FOR USE

QUARTYSEPT

READY FOR USE SOLUTION FOR FAST DISINFECTION OF SURFACES OF NON-INVASIVE MEDICAL DEVICES



CE 1639

Application

Quartysept is intended for rapid disinfection of surfaces of non-invasive medical devices and surfaces in the medical area e.g. examination tables, laboratory benches, treatment chairs, unit surfaces, instrument trays, handpieces, covers, mattresses ect.

Properties

Quartysept is a ready for use solution based on the combination of alcohols and quaternary ammonium compounds. It is aldehyde and phenol free. The preparation shows a broad spectrum of microbiological activity at short exposition time. Organic soil doesn't affect its biocide activity.

Appearance Transparent liquid
pH 5.0 - 8.0 (20°C)

Active ingredients (per 100 g):

Propan-2-ol (CAS67-63-0) – 20 g

Ethanol (CAS 64-17-5) – 10 g

Benzalkonium chloride (CAS 68424-85-1) – 0.20 g

Classification

Quartysept meets the requirements for chemical disinfectants of the MDD 93/42/CEE and is classified as a medical device - Class IIA

FOR PROFESSIONAL USE ONLY

Spectrum of activity: Bactericidal and fungicidal

Microbiological activity	Test method	Exposition
Bactericidal (<i>P.aeruginosa</i> , <i>A. S.aureus</i> , <i>E.hirae</i>)	EN 14561	1 min
Fungicidal (<i>C.albicans</i>)	EN 14562	1 min



ZHIVAS LTD



Signal word: Warning

Hazard statements

H226 - Flammable liquid and vapors. H319 Causes serious eye irritation.

H319 - Causes serious eye irritation.

H336 - May cause drowsiness and dizziness.

Precautionary statements

P102-Keep out of reach of children

P210 - Keep away from heat/sparks/open flames/hot surfaces. No smoking.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P403+P233 - Store in a well-ventilated place. Keep container tightly closed.

P337+P313 - If eye irritation persists: Get medical advice/attention.

Method of application:

Quartysept is ready for use solution. The disinfection is carried out by spraying or by wiping the surfaces of the medical devices. Exposition time is 1 minute.

Spray disinfection

Spray the surface to be treated, ensuring in particular that it is completely wetted.

Wipe disinfection

Pour the product out of the dosing bottle onto a clean cloth or directly onto the surface to be treated (approx. 25 mL product/m²). Gently wipe the surface all over with steady movements. Make sure that the surface is completely wetted during the entire exposition time.

Warnings: During the treatment of medical devices and equipment, they have to be away from sources of ignition. Do not apply on permeable and alcohol-sensitive surfaces!

Storage conditions:

30° C

In tightly closed original container in dry and well ventilated premises at t 0° C .

Treatment of waste and packaging – according to local regulations.

Do not release the concentrate in the environment – for disposal of waste, residuals and empty packaging see special instructions in safety data sheet.

Packages: (HDPE) vials of 750 ml with sprayer pump, bottles (HDPE) of 1 L and canister (HDPE) of 5 L.

Shelf life:

2 years

ZHIVAS Ltd

14 Assen Yordanov blvd., Sofia 1592, Bulgaria

Tel/fax: (+359 2) 981 78 23, (+359 2) 980 25 07

zhivas@techno-link.com;

www.zhivas.com